Haematocrit-Rotor 7500 3473

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Haematocrit-Rotor 7500 3473

D

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Haematocrit-rotor 7500 3473

UK

These instructions for use contain information about operation and cleaning of the Haematocrit-rotor.

Please read this information carefully before using the rotor.

The Haematocrit-rotor must be operated only in instrument approved for it.

Please check the manual of your centrifuge for information whether the rotor is approved for your instrument.

You will also find further information regarding the centrifuge operation in the centrifuge's manual.

Safety instructions in this manual

This symbol denotes potential hazards to persons.

This symbol denotes potential damage to the centrifuge or parts in its immediate surroundings.

General danger area

Before switching on the centrifuge please read this manual!

Technical data

Haematocrit-rotor 7500 3473

- Rotor for 24 blood capillary tubes
 - Diameter 1.2 up to 1.8 mm
 - Length maximum 75 + 1 mm
- Maximum speed 14800 min⁻¹
- Maximum RCF-value 20815 xg (at maximum speed)

Operation of the Haematocrit-rotor

Loading the Haematocrit-rotor

After filling the capillary tubes with a blood column of approximately 65 mm and closing them, insert the tubes in the rotor slots. The closed side has to face the rotor rim.

To avoid imbalance, irregular running and possible damages, the rotor always needs to be loaded symmetrically!

If the rotor is to be loaded only partially, opposing slots always need to be occupied by a capillary tube.

To close the rotor, place the rotor lid centrically on the rotor nut. By clockwise turning of the lid-handle the rotor lid is fixed.

Installing and operating the Haematocrit-rotor

After checking the correct loading of the Haematocrit-rotor, it can be installed in the centrifuge.

1. Open the centrifuge lid and make sure that the rotor chamber and the rotor are clean. Remove eventual dust, debris or sample residues. The thread and the O-Ring on the motor shaft must be in perfect condition.
2. Turn the rotor so that the notch for engaging the drive shaft points downward.
3. Place the rotor on top of the drive shaft so that the notch of the rotor is located precisely above the retaining pin. (There are two bars in the labelling on the upper side of the rotor indicating the position of the notch). These bars make the positioning easy.
4. Press the rotor gently down until it stops.
5. Grip the rotor tightly and use the provided rotor wrench to tighten up the rotor.

Do not push the rotor down using force.

If the rotor can not be tightened, you have to remove it carefully, align the notch and the retaining pin again and re-install it.

6. Place the rotor lid onto the rotor, and attend to the tight fit of the rotor.
7. Close the centrifuge lid and select speed and run time (for detailed information about selecting the parameters see the centrifuge manual). In accordance with DIN 58933, we recommend the following values to obtain reproducible results:

speed in rpm	Run time in minutes
13300	12
14800	10

8. Start the centrifuge.


Figure 1 Alignment rotor and motor

Notes on the Haematocrit determination

To obtain reproducible results, standardization is required both in the interest of patients and blood donors and to guarantee the comparability of research results.

The DIN 58933 norm describes a practical standard. The parameters for centrifugation in these notes are based on the norm.

Remove the Haematocrit-rotor

To remove the Haematocrit-rotor, follow the steps for installation in reverse order.

Measuring the Haematocrit value

The evaluation occurs after removing the rotor lid.

Reading graph 7600 0938

The reading graph comes with the rotor. To read the values, remove the capillary tubes one by one from the Haematocrit-rotor.

Place the bottom end of the capillary tubes on the zero line with the blood column and the upper end of the plasma column on the 100% line of the graph. The dividing line between erythrocytes and plasma indicates the percentage of packed cells. Be careful to place the capillary tubes at right angles to the zero line.


Figure 2 Haematocrit reading graph

for your notes

Cleaning and Care

Clean the rotor and the accessories regularly and in case of need. This is indicated both for reasons of hygiene and to prevent corrosion due to contamination sticking to the instrument and its accessories.

Clean them with mild agents of pH values ranging from 6 to 8.

Immediately after cleaning, dry the aluminum parts or put them into a warm-air dryer at a temperature not exceeding 50°C.

Disinfection

If a capillary tube containing infectious material leaks, during a centrifugation run, you have to disinfect immediately the rotor.

Infectious material could enter the centrifuge if spills or tube breakage occur.

Danger of infection may occur upon contact!

Take appropriate protective measures for personnel!

Mind the permissible filling volumes and loading limits for the tubes!

Broken glass capillaries have sharp edges and pose a risk of injury.

In case of contamination the operator has to make sure, that no further persons are jeopardized!

Contaminated parts have to be decontaminated immediately.

If required further protective measures have to be initiated.

Rotor and rotor chamber must be treated with a neutral, universal disinfectant. Best suited for this purpose are disinfectant sprays, ensuring that all rotor and accessory surfaces are covered evenly.

Please note the safety measures and handling hints when applying these substances!

You may disinfect the rotor and the accessories as described in the following section. Be sure to follow the pertinent safety procedures for handling infectious material.

1. Remove the rotor lid and the rotor from the centrifuge.
2. Remove the capillary tubes and dispose of them.
3. Treat the rotor and the rotor lid according to the instructions given for the disinfectant (soaking in liquid or spraying). You must observe the specified action times!

4. Turn the rotor head down and drain off the disinfectant. Thereafter thoroughly rinse rotor and rotor lid with water.
5. Dispose of the disinfectant according to valid regulations.
6. Aluminum rotors have to be treated with anti-corrosive protective oil subsequently (included in delivery).

Disinfection with bleaching lye

These agents contain highly aggressive hypochlorites and must not be used with aluminum rotors!

Autoclave

Rotor and rotor lid may be autoclaved at a maximum temperature of 134°C.

Maximum autoclave cycles:

- 15 minutes at 134 °C
- 20 minutes at 121 °C

Before autoclaving it, clean the rotor and rinse it with water. Remove the rotor lid and take out the capillary tubes.

Chemical admixtures in the vapour are not permitted.

Never exceed the permitted values for autoclave temperature and time.

If the rotor shows signs of wear or corrosion it may no longer be used!

If the white rubber strip shows signs of wear (e. g. cracks) it needs to be replaced.

Replace rubber strip 7500 3437

Remove the old rubber strip. Put the new rubber strip together in a ring so that it is not twisted. Put the touching ends into the rotor groove so that they are placed between two capillary slots (see figure 3). Press the rubber strip completely into the groove, avoiding creases or waves.

Figure 3 Inserting the white rubber strip

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Thermo Electron LED GmbH