## ROTANTA 460 Robotic

Ⓓ	<b>Bedienungsanleitung.....</b>	<b>5</b>
Ⓔ	<b>Operating Instructions.....</b>	<b>31</b>
Ⓕ	<b>Mode d'emploi .....</b>	<b>55</b>
Ⓘ	<b>Istruzioni per l'uso .....</b>	<b>81</b>


Fig. 1


Fig. 2 ROTANTA 460 Robotic

**EG-Konformitätserklärung  
EC Conformity Declaration  
Déclaration de conformité CE  
Dichiarazione di conformità alle norme CEE**

Andreas Hettich GmbH & Co. KG • Föhrenstraße 12 • D-78532 Tuttlingen • Germany

Das bezeichnete Gerät, inklusive Zubehör entspricht den aufgeführten EG-Richtlinien und Normen.  
The denoted device, including accessories corresponds to the listed EC guidelines and standards.  
L'appareil désigné, y compris les accessoires, correspond aux directives CE et aux normes énumérées.  
L'apparecchio designato, compresi gli accessori, è conforme alle direttive CE e alle norme citate.

Geräteart, Type of device, Type d'appareil, Tipo di apparecchio:

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## Contents

1	Use according to specification .....	33
2	Residual risks .....	33
3	Technical specifications .....	33
4	Notes on safety .....	35
5	Symbol meanings .....	36
6	Delivery checklist .....	37
7	Unpacking the centrifuge .....	38
8	Initial operation .....	38
9	Interface .....	38
10	Opening and closing the lid .....	38
10.1	Opening the lid .....	38
10.2	Closing the lid .....	39
11	Installation and removal of the rotor .....	39
12	Loading the rotor .....	39
13	Control and display elements .....	40
13.1	Control knob .....	40
13.2	Keys and setting options .....	40
14	Enter centrifugation parameters .....	41
14.1	Running time .....	41
14.2	Speed (RPM) .....	41
14.3	Relative centrifugal force (RCF) and centrifuging radius (RAD) .....	41
14.3.1	Centrifuging radius (RAD) .....	41
14.3.2	Relative centrifugal force (RCF) .....	42
14.4	Run-up and run-down parameters .....	42
14.4.1	Run-up step .....	42
14.4.2	Braking step .....	42
14.4.3	Brake switch-off speed .....	42
14.5	Temperature .....	42
15	Programming .....	42
15.1	Inputting or changing programs .....	42
15.2	Calling up programs .....	43
15.3	Write protection for programs .....	43
15.4	Program linkage .....	43
15.4.1	Activating or deactivating program linkage .....	43
15.4.2	Linking programs or changing a program linkage .....	43
15.4.3	Calling up program linkage .....	44
15.5	Automatic temporary storage .....	44
16	Centrifugation .....	44
16.1	Centrifugation with preselected time .....	44
16.2	Continuous operation .....	45
16.3	Short-time centrifugation .....	45
17	Change the settings during the centrifugation run .....	45
18	Emergency Stop .....	45

19	Acoustic signal .....	45
20	Displayed centrifugation data after switch-on.....	46
21	Set temperature unit.....	46
22	Centrifuge address.....	46
23	RCF integral and B ramp functions .....	46
24	Query of operating hours and the number of centrifugation runs .....	46
25	Querying system information.....	47
26	Immediate display of centrifugation data after switch-on.....	47
27	Cooling (only in centrifuges with cooling) .....	47
27.1	Standby-cooling .....	47
27.2	Precooling of the rotor.....	47
28	Heating (only on centrifuges with heating/cooling option) .....	47
29	Key switch.....	48
30	Setting rotor position 1 .....	48
31	Relative centrifugal force (RCF) .....	49
32	Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm <sup>3</sup> .....	49
33	Rotor Identification .....	49
34	Emergency release .....	50
35	Maintenance and servicing.....	50
35.1	Centrifuge (housing, lid and centrifuging chamber).....	50
35.1.1	Surface cleaning and care .....	50
35.1.2	Surface disinfection .....	50
35.1.3	Removal of radioactive contaminants .....	51
35.2	Rotors and Attachments .....	51
35.2.1	Cleaning and care.....	51
35.2.2	Disinfection .....	51
35.2.3	Removal of radioactive contaminants .....	51
35.2.4	Trunnions.....	51
35.2.5	Rotors and accessories with limited service lives .....	52
35.3	Autoclaving .....	52
35.4	Centrifuge containers.....	52
36	Faults .....	53
37	Acceptance of the centrifuges for repair.....	54
38	Disposal .....	54
39	Anhang / Appendix.....	107
39.1	Rotoren und Zubehör / Rotors and accessories.....	107

### 1 Use according to specification

The machine presented here is a medical product (laboratory centrifuge) according to the IVD guideline 98/79/EG. The centrifuge is used to separate substances or substance mixtures with a density of max. 1.2 kg/dm<sup>3</sup>. This also includes substances and substance mixtures of human origin. The centrifuge is only intended to be used for this purpose. A different use or application over and above this is deemed not in accordance with the specifications. The company Andreas Hettich GmbH & Co. KG undertakes no liability for damages resulting therefrom.

Belonging to the application according to specification is also the observance of all references contained in the Instruction Manual and compliance with the inspection and maintenance works.

If the centrifuge will be built into another machine or will be integrated in a system the manufacturer of the complete system is responsible for its security.

### 2 Residual risks

The machine is constructed according to the state of the art and the recognized technical safety regulations. Improper use and handling can result in dangers to life and limb of the user or third parties and impairments to the machine or to other material assets. The machine is only to be used for the specified applications and only in an impeccable technical safety condition.

Disturbances that can interfere with the safety are to be immediately rectified.

### 3 Technical specifications

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen				
Model	ROTANTA 460 Robotic				
Type	5680-70, 5680-70 (H)	5680-78, 5680-78 (H)	5680-75, 5680-75 (H)	5680-71, 5680-71 (H)	5680-74, 5680-74 (H)
Mains voltage (± 10%)	230-240 V 1~	200-220 V 1~	127 V 1~	110-120 V 1~	100 V 1~
Mains frequency	50-60 Hz				
Connected load	max. 1800 VA				
Cooling medium	R 404A				
Max. capacity	4 x 750 ml				
Allowed density	1.2 kg/dm <sup>3</sup>				
Speed (RPM)	6200				
Force (RCF)	5845				
Kinetic energy	40000 Nm				
Obligatory inspection (BGR 500)	yes				
Ambient conditions (EN / IEC 61010-1)	Indoors only				
– Set-up site	Up to 2000 m above sea level				
– Altitude	5°C to 35°C				
– Ambient temperature	Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C.				
– Humidity	II				
– Excess-voltage category (IEC 60364-4-443)	2				
– Pollution degree	I				
Device protection class	I				
Not suitable for use in explosion-endangered areas.					
EMC	EN / IEC 61326-1, Class B; FCC Class A				
– Emitted interference, Interference immunity					
Noise level (dependent on rotor)	≤ 62 dB(A)				
Dimensions	554 mm				
– Width	697 / 762* mm (* with cables connected)				
– Depth	723 mm				
– Height					
Weight	approx. 159 kg				

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen					
Model	ROTANTA 460 Robotic					
Type	5680, 5680 (H), 5680-RS232 5680-RS232 (H), 5680-10, 5680-10 (H), 5680-10-RS232, 5680-10-RS232 (H), 5680-30-A, 5680-30-A (H), 5680-30-D, 5680-30-D (H), 5680-40, 5680-40 (H)	5680-08, 5680-08 (H), 5680-08-RS232, 5680-08-RS232 (H), 5680-18, 5680-18 (H), 5680-18-RS232, 5680-18-RS232 (H), 5680-38-A, 5680-38-A (H), 5680-38-D, 5680-38-D (H), 5680-48, 5680-48 (H)	5680-05, 5680-05 (H), 5680-05-RS232, 5680-05-RS232 (H), 5680-15, 5680-15 (H), 5680-15-RS232, 5680-15-RS232 (H), 5680-35-A, 5680-35-A (H), 5680-35-D, 5680-35-D (H), 5680-45, 5680-45 (H)	5680-01, 5680-01 (H), 5680-01-RS232, 5680-01-RS232 (H), 5680-11, 5680-11 (H), 5680-11-RS232, 5680-11-RS232 (H), 5680-31-A, 5680-31-A (H), 5680-31-D, 5680-31-D (H), 5680-41, 5680-41 (H)	5680-04, 5680-04 (H), 5680-04-RS232, 5680-04-RS232 (H), 5680-14, 5680-14 (H), 5680-14-RS232, 5680-14-RS232 (H), 5680-34-A, 5680-34-A (H), 5680-34-D, 5680-34-D (H), 5680-44, 5680-44 (H)	
Mains voltage ( $\pm 10\%$ )	230-240 V 1~	200-220 V 1~	127 V 1~	110-120 V 1~	100 V 1~	
Mains frequency	50-60 Hz					
Connected load	max. 1800 VA					
Cooling medium	R 404A					
Max. capacity	4 x 750 ml					
Allowed density	1.2 kg/dm <sup>3</sup>					
Speed (RPM)	6200					
Force (RCF)	5845					
Kinetic energy	40000 Nm					
Obligatory inspection (BGR 500)	yes					
Ambient conditions (EN / IEC 61010-1)	Indoors only Up to 2000 m above sea level 5°C to 35°C Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C.					
– Set-up site						
– Altitude						
– Ambient temperature						
– Humidity						
– Excess-voltage category (IEC 60364-4-443)	II					
– Pollution degree	2					
Device protection class	I					
Not suitable for use in explosion-endangered areas.						
EMC	EN / IEC 61326-1, Class B; FCC Class A					
– Emitted interference, Interference immunity						
Noise level (dependent on rotor)	$\leq 62$ dB(A)					
Dimensions	554 mm 697 / 762* mm (* with cables connected) 684 mm					
– Width						
– Depth						
– Height						
Weight	approx. 159 kg					

#### 4 Notes on safety

**No claim of warranty will be considered by the manufacturer unless ALL instructions in this manual have been followed.**

- **The centrifuge should be installed on a good, stable base.**
- **Before using the centrifuge absolutely check the rotor for firm placement.**
- **When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.**
- **If the centrifuge will be built into another machine or will be integrated in a system the manufacturer of the complete system is responsible for its security.**
- **Rotors, suspensions and accessories that possess traces of corrosion or mechanical damage or if their term of use has expired may not be used any longer.**
- **The centrifuge may no longer be put into operation when the centrifuging chamber has safety-related damages.**
- **With swing-out rotors the trunnions must be regularly lubricated (Hettich Lubricating Grease No. 4051) in order to ensure consistent swinging out of the hangers.**

- **Before the initial operation of your centrifuge you should read and pay attention to the operating instructions. Only personnel that has read and understood the operating instructions are allowed to operate the device.**
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner. These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate. However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either heat up for at least 3 hours in the warm room before being connected to the mains, or run hot for 30 minutes in the cold room.
- Only the rotors and accessories approved by the manufacturer for this device may be used (see section "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").
- The centrifuge rotor may only be loaded in accordance with the chapter "Loading the rotor".
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm<sup>3</sup>.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with:
  - inflammable or explosive materials
  - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.  
For hazardous substances centrifuge containers with special screw caps must strictly be used. In addition to the screw cap centrifuge containers, for materials in hazard category 3 and 4 a biosafety system must be used (see the World Health Organisation's "Laboratory Biosafety Manual").  
No biosafety systems are available for this centrifuge.
- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.

- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Andreas Hettich GmbH & Co. KG company are allowed to be utilised.
- The following safety regulations apply:  
EN / IEC 61010-1 and EN / IEC 61010-2-020 as well as their national deviations.
- The safe operation and reliability of the centrifuge can only be guaranteed if:
  - the centrifuge is operated in accordance with the operating instructions,
  - the electrical installation on the site where the centrifuge is installed conforms to the demands of EN / IEC stipulations,
  - prescribed tests to BGV A1, BGR 500 are carried out by an expert.
- With centrifuges for robotic use please pay attention the notes of the key operated switch.

**5 Symbol meanings**

 Symbol on the machine:  
Attention, general hazard area.  
Before using the centrifuge implicitly read the operating instructions and pay attention to the safety relevant references!


 Symbol in this document:  
Attention, general hazard area.  
This symbol refers to safety relevant warnings and indicates possibly dangerous situations.  
The non-adherence to these warnings can lead to material damage and injury to personal.

 Symbol on the machine:  
Beware of squeezing the hands.

 Symbol on the machine and in this document:  
Beware of biohazard.

 Symbol on the machine and in this document:  
Beware of hot surface.  
Nonobservance of this warning can lead to material damage and personal injury.

 Symbol on the machine:  
Switch positions of the key-operated switch.

 Symbol on the machine:  
RS232 interface (only for centrifuges with RS232 interface).

 Symbol on the machine:  
Optical interface (only for centrifuges with optical interface).

 Symbol in this document:  
This symbol refers to important circumstances.

 Symbol on the machine and in this document:  
Symbol for the separate collection of electric and electronic devices according to the guideline 2002/96/EG (WEEE). The device belongs to Group 8 (medical devices).  
Applies in the countries of the European Union, as well as in Norway and Switzerland.

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\* The connecting cable is included in the delivery as ordered.

The rotor(s) and associated accessories are included in the delivery in the quantity ordered.

## 7 Unpacking the centrifuge

- Lift the carton upward and remove the padding.

Do not lift by the front panel.

Observe the weight of the centrifuge, refer to chapter "Technical specifications".

With a suitable number of helpers, hold the centrifuge on both sides and lift down from the pallet.

## 8 Initial operation

- According to the laboratory instrument standards EN / IEC 61010-2-020 an emergency switch to separate power supply in the event of a failure must be installed in the building electrical system. This switch has to be placed remote from the centrifuge, preferred outside of the room in which the centrifuge is installed or near by the exit of this room.
- **Position the centrifuge in a stable and level manner in a suitable place. During set-up, the required safety margin of 300 mm around the centrifuge is to be kept according to EN / IEC 61010-2-020.**

**When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.**

**If the centrifuge will be built into another machine or will be integrated in a system the manufacturer of the complete system is responsible for its security.**

- Do not place any object in front of the ventiduct.  
Keep a ventilation area of 300 mm around the ventiduct.
- Centrifuge with RS232 interface:  
Connect the RS232 centrifuge interface with an RS232 connecting cable (included in the scope of delivery for certain models) to the PC.  
Centrifuge with optical interface:  
Connect the optical centrifuge interface with a light-conducting cable (included in the scope of delivery for certain models) to the PC.
- Check whether the mains voltage tallies with the statement on the type plate.
- Connect the centrifuge with the connection cable to a standard mains socket. For connection ratings refer to Chapter "Technical specifications".
- Switch on the mains switch (switch position "I"). The LEDs in the keys are blinking.  
The following displays appear one after the other:
  1. The centrifuge model
  2. The program version
  3. The rotorcode (rotor), the maximum rotor speed (Nmax) and the preset centrifuging radius (R) of the last rotor that was recognized by rotor recognition.
  4. **⚡ OPEN ⚡ OEFFNEN**
- Open the lid.  
The centrifuge data of the last used program or program 1 will be displayed.
- Remove the transport safety device (see instruction sheet on "Moving the equipment safely").

## 9 Interface

The device is either equipped with an RS232 or an optical interface.

The RS232 interface is labeled with the  symbol.

The optical interface is labeled with the  symbol.

The centrifuge can be controlled and data queried via this interface.

The LED in the  key lights up during data communication.

## 10 Opening and closing the lid

### 10.1 Opening the lid

The lid can only be opened when the centrifuge is switched on and the rotor is at rest. If it cannot be opened under these circumstances, see the section on "Emergency release".

- Press the button . The lid unlocks via the motor and the left LED in the pushbutton  extinguishes.

## 10.2 Closing the lid


Do not put your fingers between lid and housing.  
Do not bang the lid shut.

If the left LED in the **OPEN/STOP** button flashes, press the **OPEN/STOP** button so that the motor-driven lid lock goes into the basic position (opened).

- Place the lid and lightly press down the front edge of the lid. The locking action is effected by motor. The left LED in the button **OPEN/STOP** lights up.

## 11 Installation and removal of the rotor


- Clean the motor shaft (C) and the rotor drilling (A), and lightly grease the motor shaft afterwards. Dirt particles between the motor shaft and the rotor hinder a perfect seating of the rotor and cause an irregular operation.
- Place the rotor vertically on the motor shaft. The motor shaft dog (D) has to fit in the rotor slot (B). The alignment of the groove is labelled on the rotor.
- Tighten the rotor tension nut with the supplied wrench by turning in a clockwise direction.
- Check the rotor for firm seating.
- Loosening the rotor: Loosen the tension nut by turning in a counter clockwise direction, and turning until the working point for lifting. After passing the working point for lifting the rotor is loosened from the motor shaft cone. Turn the tension nut until the rotor is able to be lifted from the motor shaft.

## 12 Loading the rotor


Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2).

- Check the rotor for firm seating.  
With swing-out rotors all rotor positions must be lined with **identical** hangers. Certain hangers are marked with the number of the rotor position. These hangers may only be used in the respective rotor position. Hangers that are marked with a set number (e.g. S001/4) may only be used in the set.
- The rotors and hangers may only be loaded symmetrically. The centrifuge containers have to be distributed evenly on all rotor positions. For authorised combinations see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

Rotor is evenly loaded

**Not permitted!**

Rotor is not evenly loaded

- On certain hangers, the weight of the maximum load or the weight of the maximum load and the maximum weight of the completely loaded hanger is specified. This weight may not be exceeded. In case of exception, see chapter "Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm<sup>3</sup>". The weight specified for the maximum loading includes the total weight of adapter, frame, centrifuging container and content.
- In containers with rubber inserts, the same number of rubber inserts must always be among the centrifuge containers.
- Always fill the centrifuge containers outside of the centrifuge.
- No liquid should be allowed to enter the centrifugal chamber during filling and swinging out of the hangers.
- The maximum filling quantity for the centrifuge containers specified by the manufacturer must not be exceeded.
- In order to maintain the weight differences within the centrifuge container as marginal as possible, a consistent fill level in the containers is to be heeded.

## 13 Control and display elements

See figure on page 2.

Fig. 2, Fig. 3: Display and control panel

### 13.1 Control knob

For setting the individual parameters.  
Turning anticlockwise reduces the value. Turning clockwise increases the value.

### 13.2 Keys and setting options

- Running time, parameters **t/hms**.  
**h**: hours. Adjustable from 1 h to 99 h, in 1 hour increments.  
**m**: minutes. Adjustable from 1 min to 59 min, in 1 minute increments.  
**s**: seconds. Adjustable from 1 s to 59 s, in 1 second increments.
- Continual running " $\infty$ "

- Revolution, parameters **rpm**.  
Adjustable from 50 rpm to a maximum rotor speed of the rotor (Nmax), in increments of 10. For maximum rotor speed, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

- Relative centrifugal force, parameter **RCF**.  
The RCF is shown in brackets  $\langle \rangle$ . The LED in the button lights up.  
The figure that is adjustable is the figure that produces a revolution speed between 50 rpm and the maximum rotor speed (Nmax). Adjustable in increments of 1.
- Centrifuging radius, parameters **RAD**.  
Adjustable from 85 mm to 330 mm in 1 millimeter increments. For centrifuging radius, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

- Run-in and run-out parameters
  - Run-in steps, parameters  $\swarrow$ .  
Step 9 = shortest run-in time, ... Step 1 = longest run-in time.
  - Braking step, parameters  $\searrow$ .  
Step 9 = shortest run-out time, ... Step 1 = long run-out time, step 0 = non-braked run-out.
  - Brake switch-off revolution, parameters **N Brake**.  
Adjustable from 50 rpm to the maximal rotor speed (Nmax) in increments of 10.  
Once the rotor speed has been reached, non-braked run-out begins.

- Temperature  
Adjustable in degrees Celsius (°C) or in degrees Fahrenheit (°F). For setting the temperature unit, see chapter "Temperatur-Einheit einstellen/Setting the temperature unit".  
Parameters **T/°C** = degrees Celsius (°C).  
Adjustable from -20°C to +40°C, in 1°C increments (for the heat/cool option, adjustable in the range from -20°C to +60°C).  
Parameters **T/°F** = degrees Fahrenheit (°F).  
Adjustable from 4°F to +104°F, in 1°F increments (for the heat/cool option, adjustable in the range from -4°F to +140°F).  
The lowest temperature that can be reached depends on the rotor (see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").

- Call up programs and program linkages; parameter **RCL**.  
Programs: Program locations 1 to 99 and PREC. program linkages: Program locations A to Z.
- Save programs and program linkages; parameter **STO**.  
99 programs can be saved. The program place PREC (PRECOOLING) is reserved for the precooling program. Program place 0 is used for temporary storage for centrifuge data from the last centrifugation run. No programs can be saved at this program place.  
It is possible to save 25 program linkages (program locations A to Z; there is no program location J). A program linkage can consist of 20 programs.
- Link programs; parameter **EDIT**.
- Start "Machine Menu" (keep the key held down for 8 seconds) and select the parameters from this menu.

- Start rotor precooling. The LED in the key is lit during the centrifugation run, as long as the rotor is revolving.  
Precooling of the rotor takes place automatically with the **PREC** (PRECOOLING) program.

- 
  - Start centrifugation run. The LED in the key is lit during the centrifugation run, as long as the rotor is revolving.
  - Short-time centrifugation.  
The centrifugation run occurs while the key is kept pressed. The LED in the key is lit during the centrifugation run, as long as the rotor is revolving.
  - Saving entries and changes.
- 
  - End centrifugation run.  
The rotor decelerates with the preset run-down parameters. The right-hand LED in the button lights up until the rotor is stationary. Once the rotor is stationary the left-hand LED flashes in the button. Pressing the button twice triggers the EMERGENCY STOP.
  - Unlock the lid.  
The left-hand LED in the button goes out.
  - Leave parameter input and the "Machine Menu".

## 14 Enter centrifugation parameters

 If no key is pressed for 8 seconds long after the selection or during the input of parameters, the previous values will be shown in the display. The input of parameter then has to be executed again.




When you enter several parameters, the  key does not have to be pressed until you have made the settings for the last parameter.

If parameters are changed, the program place number is displayed in parentheses (). This means that the centrifugation data in the display no longer corresponds to the centrifugation data from the program place that has been saved.


You can no discontinue entering parameters at any time by pressing the key . In this case the adjustments are not saved.

### 14.1 Running time

 In order to set continual running, the minutes, seconds and hours have to be set to zero. Continual running is represented in the display by the following symbol, "∞".

- Press the  key. The parameters **t/hms** is displayed. The minutes (**m**) are shown in parentheses < >, and can be changed.
- Use the adjusting knob  to set the value you want.
- Press the  key. The seconds (**s**) are shown in parentheses < > and can be changed.
- Use the adjusting knob  to set the value you want.
- Press the  key. The hours (**h**) are shown in parentheses < > and can be changed.
- Use the adjusting knob  to set the value you want.
- Press the  key, to have the settings applied to the display.

### 14.2 Speed (RPM)

- Press the  key. The parameters **RPM** is displayed.
- Use the adjusting knob  to set the value you want.
- Press the  key, to have the settings applied to the display.

### 14.3 Relative centrifugal force (RCF) and centrifuging radius (RAD)

#### 14.3.1 Centrifuging radius (RAD)

- Keep pressing the  key until the parameters **RAD** and **RCF** are displayed, e.g. **RAD = 146 -> RCF = 3695**. The LED is lit in the key.
- Use the adjusting knob  to set the centrifuging radius you want.  
By changing the centrifuging radius, the value adjusts automatically to the RCF.
- Press the  key, to have the settings applied to the display.

### 14.3.2 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) depends on the radius of the centrifuge. Before setting the RCF, the radius of the centrifuge must be set. See chapter "Centrifuge radius (RAD)".

- Keep pressing the  $\text{RCF}$  key until the parameters **RCF** and **RAD** are displayed, e.g. **RCF = 3695 RAD = 146**. The LED is lit in the key.
- Use the adjusting knob  $\odot$  to set the RCF you want.
- Press the  $\text{START}$  key, to have the settings applied to the display.

## 14.4 Run-up and run-down parameters

### 14.4.1 Run-up step

- Press the  $\text{Run-up}$  key as often as necessary until the parameter  $\text{Run-up}$  is shown.
- Use the adjusting knob  $\odot$  to set the level you want.
- Press the  $\text{START}$  key, to have the settings applied to the display.

### 14.4.2 Braking step

- Press the  $\text{Braking}$  key as often as necessary until the parameter  $\text{Braking}$  is shown.
- Use the adjusting knob  $\odot$  to set the level you want.
- Press the  $\text{START}$  key, to have the settings applied to the display.

### 14.4.3 Brake switch-off speed

- Press the  $\text{Brake switch-off}$  key as often as necessary until the parameter **N Brake** is shown.
- Use the adjusting knob  $\odot$  to set the value you want.
- Press the  $\text{START}$  key, to have the settings applied to the display.

## 14.5 Temperature

The temperature can be entered in degrees Celsius ( $^{\circ}\text{C}$ ) or in degrees Fahrenheit ( $^{\circ}\text{F}$ ). For setting the temperature unit, see chapter "Set temperature unit".

- Press the  $\text{Temp}$  key. The following parameters are displayed: **T/ $^{\circ}\text{C}$**  or **T/ $^{\circ}\text{F}$** .
- Use the adjusting knob  $\odot$  to set the value you want.
- Press the  $\text{START}$  key, to have the settings applied to the display.

## 15 Programming

If parameters are changed, the program place number is displayed in parentheses (). This means that the centrifugation data in the display no longer corresponds to the centrifugation data from the program place that has been saved.

### 15.1 Inputting or changing programs

- Set the parameters you want (see chapter "Enter centrifugation parameters").
- Press the  $\text{PROG}$  key as often as necessary until the **STO** parameter is shown.
- Use the adjusting knob  $\odot$  to set the program place you want.

If a "+" is shown after the program location, then this data is write-protected. In this case, the write protection has been removed first before saving (see the chapter on "Write protection for programs").

- Press the  $\text{START}$  key to have the settings saved to the program place you want. Briefly, **Program store ...** is displayed as confirmation.

The previous program location data will be overwritten when the new data is saved.

If "**Protected !!**" is displayed, then the data at the program location is write-protected and it will not be saved.

## 15.2 Calling up programs

- Press the **PROG** key. The parameters **RCL** is displayed.
- Use the adjusting knob  $\odot$  to set the program place you want.

If a "+" is shown after the program location, then this data is write-protected.

- Press the **START** key. Briefly, **Program recall ...** is displayed as confirmation. The centrifugation data of the selected program place is shown.

## 15.3 Write protection for programs

The programs can be protected against unintentional modification.

The write protection can be activated/deactivated as follows when the rotor is at a standstill:

- Call up the desired program (see the chapter on "Calling up programs").
- Press the **PROG** key. The parameters **RCL** is displayed.
- Press and hold the **PROG** button for eight seconds. After eight seconds, **Set Protection = 1-** (for example) appears in the display.
- Set  $\odot$  "+" or "-" with the adjusting knob.
  - + = program is write-protected,
  - = program is not write-protected.
- Press the **START** button to save the setting.

## 15.4 Program linkage

"Program linkage" can be used to link several centrifuge operations together.

A program linkage is only possible if this has been activated (parameter **Multi programs = on**; see the chapter on "Activating or deactivating program linkage").

### 15.4.1 Activating or deactivating program linkage

The program link can be activated/deactivated as follows when the rotor is at a standstill:

- Press and hold the **PROG** button for eight seconds. After eight seconds, **\*\*\* Machine Menu \*\*\*** appears in the display.
- Press the **PROG** key as often as necessary until the **->Settings** is displayed.
- Press the **START** key. The **SOUND / BELL = off/on** is displayed.
- Press the **PROG** key as often as necessary until the **Multi programs = off/on** is displayed.
- Set **off** or **on** with the adjusting knob  $\odot$ .
  - off** = program linkage deactivated,
  - on** = program linkage activated.
- Press the **START** key to save the setting. Briefly, **Store Settings ...** is displayed as confirmation, followed by **-> Settings**.
- Press the **OPEN/STOP** key once to leave the "Settings" menu or press the **OPEN/STOP** key twice to leave the "Machine Menu".

### 15.4.2 Linking programs or changing a program linkage

It is possible to save 25 program linkages (program locations A to Z; there is no program location J). A program linkage can consist of no more than 20 programs.

The speed adjustment from one program to the next one is always done in a program linkage with the start-up parameter of the next program.

There cannot be any continuous operation programs in a program linkage.

No centrifuge parameters may be modified in a program linkage. A parameter modification is only possible in the individual programs.

The **TIME** button can be used during the centrifuge operation to call up the total run time of the program linkage

(e. g.  $\Sigma=00:05:30$ ) and the run time of the currently running program (e.g. **t B.02=00:01:00**).

1. Press the **PROG** button repeatedly until the parameter **EDIT A...Z** is displayed.
2. Use the  $\odot$  adjusting knob to set the desired program location where the program linkage is to be saved.
3. Press the **START** button. The program location of the program linkage and the first program of the program linkage will be shown (e.g. **EDIT B.01 = 01**).
4. Use the  $\odot$  adjusting knob to set the first program of the program linkage.
5. Press the **PROG** button. The next program of the program linkage will be shown (e.g. **EDIT B.02 = END**).

6. Use the  adjusting knob to set the next program of the program linkage.
7. Press the  button. The next program of the program linkage will be shown (e.g. **EDIT B.03 = END**).
8. Repeat steps 6 and 7 until all programs have been set.
9. Set  **END** with the adjusting knob (turn adjusting knob anti-clockwise).

For program linkages which consist of 20 programs, **END** cannot be set after the 20th program.

10. Press the  button. **STO B**, for example, is displayed.
11. Press the  button to save the program linkage.  
As confirmation, **Multi program store ..** will be shown briefly.

### 15.4.3 Calling up program linkage

- Press the  button repeatedly until the parameter **RCL A...Z** is displayed.
- Use the adjusting knob  to set the program place you want.  
Press the  key. Briefly, **Multi program recall ...** is displayed as confirmation.
- The centrifugation data of the first program of the program linkage will be shown.

### 15.5 Automatic temporary storage

The program place 0 serves as temporary storage for centrifugation data of the last centrifugation run that took place.

No programs can be saved at this program place.

After every start to a centrifugation run, the centrifugation data that is used for the run is automatically saved at program place "0" and can be accessed there.

## 16 Centrifugation

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

If the centrifuge will be built into another machine or will be integrated in a system the manufacturer of the complete system is responsible for its security.

If the permissible weight difference within the rotor loading has been exceeded, the drive shuts down during the start-up, and **IMBALANCE** is displayed.

If the speed in the selected program is higher than the maximum speed of the rotor ( $N_{max}$ ), it will not be possible to start a centrifugation run. The message **N > ROTOR MAX** will be displayed (see Chapter "Faults").

If **N > ROTOR MAX in Prog** : e.g. **5**, **Runtime 00:00 in Prog** : e.g. **5** or **Empty Program** is shown during program linkage, no centrifuge operation can be started (see the chapter on "Faults").

A centrifugation run can be stopped at any time by pushing the key  .

All parameters can be selected and altered during the centrifugation run (see Chapter "Enter centrifugation parameters").

It is possible to switch over between RPM and RCF display at any time using the  and  keys. Switching is not possible if program linkages are being used. When working with the RCF display, it is necessary to enter the centrifuging radius.

If  **OPEN**  **OEFFNEN** is displayed, a further operation of the centrifuge is only possible after opening the lid once.

If the rotor has been changed, no centrifugation run takes place and this display is shown, e.g. **Rotor 4 Nmax= 4500 R=184 mm** (see chapter "Rotor Identification").

Operation errors and malfunctions will be shown (see Chapter "Faults").

- Switch on the mains supply switch. Switch position **I**.
- Load the rotor and close the centrifuge lid.

### 16.1 Centrifugation with preselected time

- Set a run time, call up a program with preset time or a program linkage (see the chapters on "Enter centrifugation parameters", "Calling up programs" or "Program linkage").
- Press the  key. The LED in the  key blinks until the rotor has been imported, it is subsequently lit.
- Once the time has run up or if the centrifugation run is canceled by pressing the  key, run-down takes place with the selected run-down parameter. The run-down parameter is displayed e.g. **9.9**. The right-hand LED in the  key is lit. Once the rotor has come to a standstill, the LED in the  key and the right LED in the  key switch off and the left LED in the  key lights up.

During the centrifugation run, the speed of the rotor or the RCF value resulting from this, the sample temperature and the remaining time will be displayed.

## 16.2 Continuous operation

- Set the minutes, seconds and hours to "0" or start a continual run program (see chapter "Enter centrifugation parameter" or "Calling up programs").
- Press the **START** key. The LED in the **START** key blinks until the rotor has been imported, it is subsequently lit. The time count starts from 00:00.
- Press the **OPEN/STOP** key to end the centrifugation run. Run-down is carried out with the selected run-down parameters. The run-down parameter is displayed e.g. 9. The right-hand LED in the **OPEN/STOP** key is lit. Once the rotor has come to a standstill, the LED in the **START** key and the right LED in the **OPEN/STOP** key switch off and the left LED in the **OPEN/STOP** key lights up.

During the centrifugation run, the speed of the rotor or the RCF value resulting from this, the sample temperature and the time run will be displayed.

## 16.3 Short-time centrifugation

A short-time centrifugation is not possible if program linkages are being used.

- Keep the **START** key pressed. The LED in the **START** key blinks until the rotor has been imported, it is subsequently lit. The time count starts from 00:00.
- Release the **START** key again to end the centrifugation run. Run-down is carried out with the selected run-down parameters. The run-down parameter is displayed e.g. 9. The right-hand LED in the **OPEN/STOP** key is lit. Once the rotor has come to a standstill, the LED in the **START** key and the right LED in the **OPEN/STOP** key switch off and the left LED in the **OPEN/STOP** key lights up.

During the centrifugation run, the speed of the rotor or the RCF value resulting from this, the sample temperature and the time run will be displayed.

## 17 Change the settings during the centrifugation run

It is not possible to modify settings during centrifuge operation if program linkages are being used.

The running time, the speed, the relative centrifugal force (RCF), the centrifuging radius, the run-up and run-down parameters as well as the temperature can be changed while the centrifuge is running.

- Select the parameter you want and use the  $\odot$  adjusting knob to change the value (see chapter "Enter centrifugation parameter").
- Press the **START** key. The changed setting is saved at program place "0" (see chapter "Automatic temporary storage"). The program place number is displayed in parentheses (0). The original program is not overwritten.

## 18 Emergency Stop

- Press the **OPEN/STOP** key twice.

During a not-stop, run-out is carried out with braking level 9 (shortest run-out). Braking level 9 is displayed. If braking level 0 has been preselected, run-out is carried out with braking level 9d. With braking level 9d, the run-out time is longer, technically speaking, than braking level 9.

## 19 Acoustic signal

The acoustic signal sounds:

- Upon the appearance of a disturbance in 2 second intervals.
- After completion of a centrifugation run and rotor standstill in 30 second intervals.

The acoustic signal is stopped by opening the lid or pressing any key.

The signal after completion of the centrifugation run can be activated or deactivated in the following manner, if the rotor is at standstill:

- Keep the **PROG** key held down for 8 seconds.  
After 8 seconds, **\*\*\* Machine Menu \*\*\*** appears in the display.
- Press the **PROG** key as often as necessary until the **->Settings** is displayed.
- Press the **START** key. The **SOUND / BELL = off/on** is displayed.
- Use the adjusting knob  $\odot$  to set **off** or **on**.
- Press the **START** key to save the setting.  
Briefly, **Store Settings ...** is displayed as confirmation, followed by **-> Settings**.
- Press the **OPEN/STOP** key once to leave the "Settings" menu or press the **OPEN/STOP** key twice to leave the "Machine Menu".

## 20 Displayed centrifugation data after switch-on.

After switch-on the centrifugation data from program 1, or from the last program that was used, is displayed. If the rotor is at standstill, this can be set as follows:

- Keep the **PROG** key held down for 8 seconds.  
After 8 seconds, **Machine Menu** appears in the display.
- Press the **PROG** key as often as necessary until the **->Settings** is displayed.
- Press the **START** key. The **SOUND / BELL = off/on** is displayed.
- Press the **PROG** key as often as necessary until **Start program = Last/First** is displayed.
- Use the adjusting knob  $\odot$  to set **Last** or **First**.  
Last = the program used last, First = program 1.
- Press the **START** key to save the setting.  
Briefly, **Store Settings ...** is displayed as confirmation, followed by **-> Settings**.
- Press the **OPEN / STOP** key once to leave the "Settings" menu or press the **OPEN / STOP** key twice to leave the "Machine Menu".

## 21 Set temperature unit

The temperature can be entered in degrees Celsius (°C) or in degrees Fahrenheit (°F). When the rotor is at standstill, the temperature unit must be set as follows:

- Keep the **PROG** key held down for 8 seconds.  
After 8 seconds, **Machine Menu** appears in the display.
- Press the **PROG** key as often as necessary until the **->Settings** is displayed.
- Press the **START** key. The **SOUND / BELL = off/on** is displayed.
- Press the **PROG** key as often as necessary until **Temp Unit = Celsius/Fahrenheit** is displayed.
- Use the adjusting knob  $\odot$  to set **Celsius** or **Fahrenheit**.
- Press the **START** key to save the setting.  
Briefly, **Store Settings ...** is displayed as confirmation, followed by **-> Settings**.
- Press the **OPEN / STOP** key once to leave the "Settings" menu or press the **OPEN / STOP** key twice to leave the "Machine Menu".

## 22 Centrifuge address

 Ex works, the address is ] = 29. Address is set.

## 23 RCF integral and B ramp functions

 These functions are not activated as standard.

## 24 Query of operating hours and the number of centrifugation runs

The operating hours are divided up into internal and external operating hours.

Internal operating hours: Total time the device was switched on.




External operating hours: Total time of previous centrifugation runs.

The operating hours and the number of centrifugation runs can be queried as follows when the rotor is at standstill:

- Keep the **PROG** key held down for 8 seconds.  
After 8 seconds, **Machine Menu** appears in the display.
- Press the **PROG** key as often as necessary until the **-> Operating Time** is displayed.
- Press the **START** key. The external operating hours are displayed, e.g. **OP Time ext = 0h25m**.
- Press the **PROG** key. The internal operating hours are displayed, e.g. **OP Time int = 1h36m**.
- Press the **PROG** key. The number of centrifugation runs is displayed, e.g. **Number of Starts = 10**.
- Press the **OPEN / STOP** key to leave the menu, "Operating Time" or press the key **OPEN / STOP** three times to leave the "Machine Menu".

## 25 Querying system information

Rotor information and the program version of the centrifuge as well as the frequency converter can be queried as follows when the rotor is at a standstill:

- Keep the  key held down for 8 seconds.  
After 8 seconds, **\*\*\* Machine Menu \*\*\*** appears in the display.
- Press the  key as often as necessary until the **->Info** is displayed.
- Press the  key. The rotorcode (rotor), the maximum rotor speed (Nmax) and the preset centrifuging radius (R) of the last rotor that was recognized by the rotor recognition, e.g. **Rotor 4\*: Nmax = 4500 R=184**.  
The last recognized rotor is marked with a star (\*). Use the adjusting knob  to have the information of the rotors permitted in the centrifuge displayed.
- Press the  key. The program version of the centrifuge is displayed, e.g. **SW version = 01.00**.
- Press the  key. The program version of the frequency converter is displayed e.g. **FC-SW-Version = 4**.
- Press the  key twice to leave the "Info" menu or press the key  three times to leave the "Machine Menu".

## 26 Immediate display of centrifugation data after switch-on.

- Switch on the power supply. (Switch setting I).
- The first time the display changes optically, (inverse display) press any key and keep it held down. The centrifugation data is displayed immediately.

## 27 Cooling (only in centrifuges with cooling)

The temperature set-point can be adjusted from -20°C to +40°C / -4°F to +104°F. On centrifuges with a heating/cooling option, the temperature set-point is settable from -20°C to +90°C / -4°F to +194°F. The lowest obtainable temperature is dependent on the rotor (see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").

### 27.1 Standby-cooling

When the rotor is at a standstill and when the lid is closed, the centrifugal chamber is cooled down to the preselected temperature if this is lower than 20°C / 68°F.

During the standby cooling, the preset temperature is displayed.






### 27.2 Precooling of the rotor

For quick precooling of the unloaded rotor and the accessories, it is a good idea to do a centrifugation run with the continuous operation setting and a speed of about 20% of the maximum rotor speed.

The precooling is done automatically with the **PREC (PRECOOLING)** program.

Precooling is not possible if program linkages are being used.

- Press the  key. The LED in the  key blinks until the rotor has been imported, it is subsequently lit. The time count starts from 00:00.
- Press the  key to end the centrifugation run. Run- down is carried out with the selected run- down parameters. The run-down parameter is displayed e.g. 9. The right-hand LED in the  key is lit. Once the rotor has come to a standstill, the LED in the  key switches off and **OPEN OEFFNEN** is displayed. The right-hand LED in the  key also switches off, the left LED in the  key starts blinking and continues blinking until the lid is opened.

During the centrifugation run, the speed of the rotor, or the resulting RCF value, the sample temperature and the remaining time or expired time are displayed.

## 28 Heating (only on centrifuges with heating/cooling option)

During operation of the centrifuge the centrifuge chamber can be heated as required to the preselected temperature. The heating is switched off when the rotor is stationary.


Risk of burning! The surface temperature of the heating element in the centrifugal chamber of the centrifuge can reach 500°C / 932°F. Do not touch the heating element.

## 29 Key switch

Only authorized skilled personnel may operate the device in the "TEACH" key setting.  
The key is to be kept in a safe place so that it is protected from unauthorized access.

The following program locks and functions can be set with the key switch when the rotor is at a standstill:

Key setting	Program locking / function
0	The centrifuge can be operated via the control panel. The centrifuge cannot be controlled via the interface.
LOCK 2	<b>LOCK 2</b> is displayed in the "↘" field. The centrifuge can be controlled via the interface. The centrifuge cannot be operated via the control panel.
TEACH	Display e.g. <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> <p><b>Teach</b></p>  </div> <div style="text-align: center;"> <p><b>Open</b></p>  </div> <div style="text-align: center;"> <p><b>=Close</b></p>  </div> </div> <p>= : Position of the hatch   : Open the hatch (keep key pressed)   : Close the hatch (keep key pressed)   : Set rotor position 1. See the chapter "Setting rotor position 1".</p>

## 30 Setting rotor position 1



Only authorized skilled personnel may set rotor position 1.

Every time the rotor is installed, rotor position 1 must be reset afterwards. Rotor position 1 is the position where place 1 of the rotor must be located for loading and unloading. Rotor position 1 can be set as follows when the rotor is at a standstill:

- Look through the window in the lid to make sure that the rotor is at a standstill.
- Turn the key to the "TEACH" setting. For example, **Teach Open =Close** might be displayed.
- Keep the  key pressed until the hatch is completely open. **Teach =Open Close** is displayed.
- Press the  key. The saved rotor position 1 is displayed, e.g. **Teach Pos 1 : 1812 <Start=ok>**.
- Turn the rotor slowly until you hear an acoustic signal. Afterwards, continue to turn it until place 1 of the rotor is in the desired loading/unloading position. Hold the rotor in this position.
- Press the  button to save the setting.  
As a confirmation, **Store Pos 1 ...** is briefly displayed, and then e.g. **Goto Pos 1/24 2826 0**.

If "**no zero pulse !**" is briefly displayed, the rotor was turned to the loading/unloading position before the acoustic signal sounded. In this case, repeat the rotor position 1 setting.

- To test the positioning, set the desired rotor place with the rotary knob , and then press the  key. The rotor continues to turn until the set rotor place is in the loading/unloading position.
- Press the  key to exit the menu. **Teach =Open Close** is displayed.
- Keep the  key pressed until the hatch is completely closed. **Teach Open =Close** is displayed.
- Turn the key back to the "LOCK 2" position.

### 31 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left( \frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge.

For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

### 32 Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm<sup>3</sup>

When centrifuging with maximum revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm<sup>3</sup>.

The speed must be reduced for materials or mixtures of materials with a higher density.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed } (n_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density [kg/dm}^3]}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, density 1.6 kg/dm<sup>3</sup>

$$n_{\text{red}} = \sqrt{\frac{1.2 \text{ kg/dm}^3}{1.6 \text{ kg/dm}^3}} \times 4000 \text{ RPM} = 3464 \text{ RPM}$$

In the exceptional case that the maximum loading indicated on the hanger is exceeded, the speed must also be reduced.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed } (n_{\text{red}}) = \sqrt{\frac{\text{maximum load [g]}}{\text{actual load [g]}}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, maximum load 300 g, actual load 350 g

$$n_{\text{red}} = \sqrt{\frac{300 \text{ g}}{350 \text{ g}}} \times 4000 \text{ RPM} = 3703 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

### 33 Rotor Identification

After every start of a centrifugation run the rotor utilised is identified.

After changing a rotor, the engine switches off and the rotorcode (rotor), the maximum rotor speed of the rotor (Nmax) and the preset centrifuging radius (R) of the rotor being used are displayed, e.g. **rotor 4 Nmax= 4500 R=184 mm**.

- Press the **OPEN/STOP** key to open the lid, or press the **START** key to start the centrifugation run. For centrifuges with refrigeration, precooling of the rotor can also be started by pressing the **☒** key.

If, following a rotor change, the maximum speed of the rotor is less than the set speed, the speed is limited to the maximum speed of the rotor. In this case, the program place number is displayed in parentheses ().

### 34 Emergency release

During a power failure the lid cannot be unlocked by motor. An emergency release has to be executed by hand.

For emergency release disconnect the centrifuge from the mains.  
Open the lid only during rotor standstill.

See figure on page 2.

- Switch off the mains switch (switch position "0").
- Look through the window in the lid to be sure that the rotor has come to a standstill.
- Insert the hexagonal wrench key into the bore hole (Fig. 1, A) and carefully rotate by half a turn in clockwise direction until the lid can be opened.
- Pull the hexagon socket head wrench out of the drilling again.
- If the left LED in the **OPEN/STOP** button flashes after the centrifuge is switched on again, press the **OPEN/STOP** button so that the motor-driven lid lock goes into the basic position (opened) again.

### 35 Maintenance and servicing

The device can be contaminated.

Pull the mains plug before cleaning.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Centrifuges, rotors and accessories must not be cleaned in rinsing machines.
- They may only be cleaned by hand and disinfected with liquids.
- The water temperature must be between 20 – 25°C.
- Only detergents/disinfectants may be used which:
  - have a pH between 5 - 8
  - do not contain caustic alkalis, peroxides, chlorine compounds, acids and alkaline solutions
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

#### 35.1 Centrifuge (housing, lid and centrifuging chamber)

##### 35.1.1 Surface cleaning and care

- Clean the centrifuge housing and the centrifuging chamber regularly, using soap or a mild detergent and a damp cloth if required. For one thing, this services purposes of hygiene, and it also prevents corrosion through adhering impurities.
- Ingredients of suitable detergents:  
soap, anionic tensides, non-ionic tensides.
- After using detergents, remove the detergent residue by wiping with a damp cloth.
- The surfaces must be dried immediately after cleaning.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- Lightly rub the rubber seal of the centrifuge chamber with talcum powder or a rubber care product after each cleaning.
- The centrifuging chamber is to be checked for damage once a year.

If damage is found which is relevant to safety, the centrifuge may no longer be put into operation. In this case, notify Customer Service.

##### 35.1.2 Surface disinfection

- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately.
- Ingredients of suitable disinfectants:  
ethanol, n-propanol, isopropyl alcohol, glutardialdehyde, quaternary ammonium compounds.
- After using disinfectants, remove the disinfectant residue by wiping with a damp cloth.
- The surfaces must be dried immediately after disinfecting.

### 35.1.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for removing radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants: anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove the agent residue by wiping with a damp cloth.
- The surfaces must be dried directly after removing the radioactive contaminants.

## 35.2 Rotors and Attachments

### 35.2.1 Cleaning and care

- In order to prevent corrosion and material changes, rotors and accessories must be cleaned regularly with soap or a mild detergent and a damp cloth. Cleaning is recommended at least once a week. Contaminants must be removed immediately.
- Ingredients of suitable detergents: soap, anionic tensides, non-ionic tensides.
- After using detergents, remove detergent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after cleaning.
- Angle rotors, container and hanger made of aluminium are to be lightly greased after drying using acid-free grease, e.g. vaseline.
- In order to prevent corrosion as a result of moisture between the rotor and the motor shaft, the rotor should be disassembled and cleaned at least once a month, and the motor shaft should be lightly greased.
- The rotors and the accessory parts have to be checked on a monthly basis for wear and corrosion. For swing-out rotors, it is important to check the area of the lifting lugs, for hangers, the grooves and the base should be checked for cracks.

Rotors and attachments may no longer be utilised upon indication of wear and tear or corrosion.

- Check the firm seating of the rotor on a weekly basis.

### 35.2.2 Disinfection

- If infectious material should get on the rotors or accessories, they must be appropriately disinfected.
- Ingredients of suitable disinfectants: glutaraldehyde, propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove disinfectant residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after disinfection.

### 35.2.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for the removal of radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants: anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove agent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after removing the radioactive contaminants.

### 35.2.4 Trunnions

With swing-out rotors the trunnions must be regularly lubricated (Hettich Lubricating Grease No. 4051) in order to ensure consistent swinging out of the hangers.

### 35.2.5 Rotors and accessories with limited service lives

The use of certain rotors, hangers and accessory parts is limited by time.

These are marked with the maximum permitted number of operating cycles or with an expiration date and the maximum permitted number of operating cycles or just with the expiration date; e.g.:

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For safety reasons, rotors, hangers and accessory parts may no longer be used if either the indicated maximum number of operating cycles or the indicated expiration date has been reached.

The number of centrifuge operations can be queried; see the chapter on "Querying the operating hours and the number of centrifuge operations".

### 35.3 Autoclaving

Swing-out rotors, angle rotors made of aluminium, suspension made of metal, lids with bio-containment as well as stands and reductions can be autoclaved at 121° C / 250°F (20 mins).

Otherwise you must ask the manufacturer.

No statement can be made about the degree of sterility.

Autoclaving accelerates the ageing process of plastics. In addition, autoclaving may discolour plastics.

### 35.4 Centrifuge containers

- With leakiness or after the breakage of centrifuging containers broken container parts and leaked centrifugation material are to be completely removed.
- The rubber inserts as well as the plastic sleeves of the rotors are to be replaced after a glass breakage.

Remaining glass splinters cause further glass breakage!

- If this concerns infectious material, a disinfection process is to be executed immediately.

### 36 Faults

If the fault cannot be eliminated with the help of the fault table, please inform Customer Service.

Please specify the type of centrifuge and the serial number. Both numbers can be found on the name plate of the centrifuge.

Perform a MAINS RESET:

- Switch off the mains switch (switch position "0").
- Wait at least 10 seconds and then switch on the mains switch again (switch position "1").

Message / fault		Cause	Remedy
No display		No voltage. Overvoltage protection tripped out.	- Check supply voltage. - Mains switch ON.
TACHO - ERROR	1, 2, 96	Faulty speedometer. Motor, electronics defective.	- Open the cover. - Switch off the mains switch (switch position "0"). - Wait at least 10 seconds. - Turn the rotor vigorously by hand. - Switch on the mains switch again (switch position "1"). The rotor must turn during switch-on.
IMBALANCE	3 *	The rotor is unevenly loaded.	- Open lid. - Check the loading of the rotor, see chapter "Loading the rotor". - Repeat the centrifugation run.
CONTROL - ERROR	4.1 - 4.5, 6	Error in lid locking or lid closure.	- Perform a MAINS RESET.
CONTROL - ERROR	4.6 - 4.8	Error / defect control of the hatch	
N > MAX	5	Rotation too fast	
N < MIN	13	Rotation too slow	
ROTORCODE	10.1, 10.3	Incorrect rotor coding	
MAINS INTERRUPT	11 *	Power failure during the centrifugation run. (The centrifugation run was not finished.)	- Open lid. - Push <b>START</b> button. - Repeat the centrifugation run if necessary.
VERSION-ERROR	12	Mismatch between electronic components Error / defect electronics	- Perform a MAINS RESET.
CONTROL-ERROR	22, 25.1-25.4	Error / defect electronics	
SER I/O - ERROR	31, 34, 36	Error / defect electronics	
POS-ERROR	40.1-47	Error / defective positioning or hatch drive	
° C * - ERROR	51, 53 - 55, 97, 98	Error / defect electronics	
° C * - ERROR	52	Temperature is too high in overspeed room. Error / defect electronics	
LOCK_ERROR	57.0, 57.1	Error / defective program locking	
FU / CCI - ERROR	60, 61.2-61.20, 61.128-61,131, 62	Error / defect electronics / motor	
FU / CCI - ERROR	61.1	Supply voltage is too low Error / defect electronics / motor	

\* Error number do not appear in the display.

Message / fault		Cause	Remedy
SENSOR-ERROR	90	Error / defect electronics	– Perform a MAINS RESET.
SENSOR-ERROR	91 - 93	Error / defect unbalance sensor	
NO ROTOR OR ROTORCODE ERROR	---	No rotor installed. Faulty speedometer.	– Open lid. – Install rotor.
N > ROTOR MAX	---	Speed in the selected program greater than the maximum speed of the rotor.	– Check the set speed. Reduce the set speed
N > ROTOR MAX in Prog : e.g. 3	---	There is a program in the displayed program location whose speed is higher than the maximum speed of the rotor.	
Runtime 00:00 in Prog : e.g. 3	---	There is a continuous operation program in the displayed program location.	– Replace the continuous operation program in the program linkage with a program with preset time.
Empty Program	---	There is no program linkage saved on the displayed program location.	– Call up a program linkage.
FC INIT ERROR	---	Error / defect electronics	– Perform a MAINS RESET.
FC VERSION ERROR		Error / defect electronics	
FATAL EEPROM ERROR	1 - 5	Error / defect electronics	

### 37 Acceptance of the centrifuges for repair

If the centrifuge is returned to the manufacturer for repair, it must be decontaminated and cleaned to protect persons, environment and material.

We reserve the right to accept contaminated centrifuges.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

### 38 Disposal

Before disposal, the device must be decontaminated and cleaned to protect people, the environment and property.

When you are disposing of the device, the respective statutory rules must be observed.

Pursuant to guideline 2002/96/EC (WEEE), all devices supplied after August 13, 2005 may not be disposed as part of domestic waste. The device belongs to group 8 (medical devices) and is categorized in the business-to-business field.

The icon of the crossed-out trash can shows that the device may not be disposed as part of domestic waste.

The waste disposal guidelines of the individual EC countries might vary. If necessary, contact your supplier.

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All rights reserved. No part of this publication may be reproduced without the prior written permission of the copyright owner.