



Calibration Laboratory and Location of Calibration

Certificate number: 020-191023

DJB Labcare Ltd 12 Howard Way Cromwell Business Centre Newport Pagnell Bucks MK16 9QS

> Tel: 01908 612598 Fax: 01908 217974

E-mail: service@djblabcare.co.uk Website: www.djblabcare.co.uk

Certificate of Calibration

Equipment: Fastime 3 Stopwatch

Serial number: 020

The results refer only to the equipment listed above

Date of Issue: 19/10/2023 Date of Calibration: 19/10/2023

Ambient Temperature: 22.5°C

Results:

Cot Value	Management Value
Set Value	Measured Value
05:00.00	05:00.04
	-
10:00.00	10:00.05
60:00.00	60:00.07

Method:

A time interval standard is used and compared to the DUT's display, carried out at 5,10 and 60 min intervals.

Uncertainty: 0.48s

Uncertainty is based on a standard uncertainty multiplied by a coverage factor k=2, providing a level of confidence of approximately 95%. The uncertainty evaluation has been carried out in accordance with UKAS requirements.

Test Equipment Details: BT Speaking Clock

Approved by: C Bradnam

END OF RESULTS