

F19-6x200 mL

Fixed Angle Carbon Fiber Rotor

6x200 mL




for
**Sorvall and Beckman
Ultra Centrifuges**

Features

- Lifetime warranty
- Safer than metallic rotors
- Corrosion free
- Cold room compatible
- Easy to handle
- Fully autoclavable
- Reduced wear on centrifuge drive
- Rapid acceleration and deceleration
- Highly resistant to radiolabel contamination

Applications

- General-purpose rotor for extra large volume
- Processing of particulate materials and debris in solution and bacteria from cell cultures
- Also subcellular organelles from tissue homogenates



Specifications

- Max. Speed: 19,000 rpm
- Max. g-force: 54,988 x g
- Max. radius (r_{max}): 136 mm
- Min. radius (r_{min}): 38 mm
- Tube Angle: 23°
- k factor: 894
- Rotor Capacity: 6 x200 mL (1200 mL)
- Max Tube Size: $\phi 62 \times 120$ mm (2.5 x 5.5 in.)
- Rotor Mass: 20 lbs
- Liquid Containment: No

Unique Features:

- Highest capacity lightweight ultra speed rotor
- General-purpose pelleting

For Sorvall Ultra Centrifuge

For Beckman Ultra Centrifuge

FIBERLite Rotor Model	F19L-6x200 mL	F19L-6x200 mL
FIBERLite Part No.	096-067030	096-067030
Will Spin In	OTD/Discovery/Ultra Pro 80	L Series
Will Replace	A-621	Type 19
Rotor Accessories and Replacement Parts		
Lid Assembly	099-067030	099-067030
O-Ring Refresher Kit	021-067030	021-067030

Ultra

Recommended Bottles with Cap Assemblies

Tube Style Material	Part No.
Polyallomer (Set of 2)	010-0372