

Cyto Rotor: New and only from Hettich

Reliable protection against infectious aerosols



Cyto Rotor





View into the centrifuging chamber showing the 1515-A rotor with 12 cyto inserts

GOOD REASONS TO CHOOSE HETTICH PRODUCTS.

LARGEST SELECTION OF CYTO ACCESSORIES

Hettich has a bigger line of cytology accessories than any other manufacturer. All of the common methods for the production of cyto preparations can now be carried out in Hettich laboratory centrifuges.

SUBSTANTIAL COST SAVINGS

Cytology laboratories that have used alternative products to date (such as Shandon® centrifuges) can now use just a single Hettich centrifuge. There is no need to acquire a centrifuge solely for cytology purposes. Users who are already Hettich customers can use the new rotors in their existing ROTOFIX 32 A and UNIVERSAL 320 models or their predecessors.

MAXIMUM FLEXIBILITY

The 12-place cyto are easy to insert and remove without the need for any tools.

THE NEW ACCESSORIES AT A GLANCE



12-place rotor 1515-A lid with biocontainment*



Clip 1524 (set of 4)



MAXIMUM SAFETY

Hettich is the first and only manufacturer to offer a cyto rotor with biocontainment. This means optimum protection against infectious aerosols.

* Tested by the TÜV in conformity with DIN EN 61010, section 2-020



Cyto chamber 1534 (500 x)

Cyto chamber 1531

(50 x)with white filter card for sample quantities in the range 200 - 500 µl



Cyto chamber 1532

(50 x)with brown filter card for very small sample quantities up to 200 µl



Cyto chamber 1538

(12 x)

(reusable) for sample quantities up to 500 µl

Filter card 1539

(200 x)for 1538



Cyto chamber 1533 (500 x)

Cyto chamber 1535

(100 x)

for volumes up to 6 ml (e.g., urine samples)



Cyto chamber 1536

(12 x)

(reusable) for volumes up to 6 ml

(e.g., urine samples)

Seal 1537

(100 x)for 1536





